

1. Our PE ponchos specifications:

MOQ	5,000 pieces
Material	PE
Size	90 * 120, 102 * 128cm (customize)
Weight	11kg / piece
Color	black
Thickness	60g PE
Features	Emossed
Packaging	OPP, 250pcs / 51 * 30 * 23 cm / piece
Lead time	11 days

2. Our PE ponchos specifications:

Our PE ponchos specifications are as follows:

Material: PE, Thickness: 60g, Weight: 11kg / piece, Color: black, Features: Emossed

3. Our PE ponchos specifications:

1. Our PE ponchos specifications are as follows: Material: PE, Thickness: 60g, Weight: 11kg / piece, Color: black, Features: Emossed

2. Our PE ponchos specifications are as follows: Material: PE, Thickness: 60g, Weight: 11kg / piece, Color: black, Features: Emossed

necessary specifications: PE / 60g, 11kg / piece, black, Emossed, OPP, 250pcs / 51 * 30 * 23 cm / piece & amp; 11 days lead time ! Foled delivered.

4. Our PE ponchos specifications:



orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INTERNATIONAL

