

Technical Specifications:

- Material: Non-woven fabric
- Weight: 30-60 gsm
- Closure: Velcro
- Usage: Laboratory, clinical, laboratory
- Size: S / M / L / XL / XXL / XXXL
- Color: White
- Quantity: 1pcs / polybag, 50pcs / Carton
- Material Composition: 80% Polypropylene
- Certification: CE, ISO13485

Size Chart:

CM	S	M	L	XL	XXL	XXXL
Length	115	115	120	125	130	135
Width	127	137	140	145	150	155

Standards:

- A.standard
- B.standard
- C.half
- D.full

Features:

- 1.SMS: Soft, comfortable, comfortable
2. Breathability
3. Provide
4. Disposable
- 5.SMS .Disposable
- 6.SMS: "SMS" (Soft, comfortable, comfortable)

