

[headwear cross kf94 mask adults multicolored](#)

All masks are made of the highest quality materials, providing excellent protection and superior comfort:

- Non-linting materials are smooth and comfortable to eliminate skin irritation
- Highly breathable construction layers remain cool, even during long cases
- Super soft earloop makes comfortable wearing for long time

[headwear cross kf94 mask adults multicolored](#)

Material	Layer 1: PP Non Woven (Polypropylene) Layer 2: Melted Filter Layer 3: Melted Filter Layer 4: PP Non Woven
Structure	4ply
Size (IZE)	20*8 CM
Weight (W)	30GSM + 25GSM + 25GSM+25GSM, C OULD Be
Standard	CE, ISO13485, FDA
Moq:	20000pcs
Color	Multicolored
Quantity	5 0 PCS/B OX, 200 0 S/CTN
Lead Time	1 Week
Delivery	7-30 Days
Payment	T/T, L/C, D/P, D/A, Western Union, Money Order, etc.

[Cross KF94 MASK multicolored](#)



orient Honest
MEDICAL INTERNATIONAL



orient Honest
MEDICAL INTERNATIONAL





 orient Honest
MEDICAL INTERNATIONAL

