

Technical Specifications:

- 1. Material: PP50g
- 2. Color: White
- 3. Weight: 50g / 1.76oz
- 4. Size: S-XXXL
- 5. Quantity: 5 pieces / polybag, 50PCS / 1.76lb
- 6.1 * 20 "FCL: 50000pcs (1.76lb)
- 7. Mini: 10 000PCS
- 8. Certifications: CE, FDA, ISO13485

Product Features:

- 1. Hand knitted
- 2. Non-breathability
- 3. Against
- 4. Usually

LAB CHART

CM	S	M	L	XL	XXL	XXXL	XXXXL	XXXXXL
100	103	107	110	112	113	120	120	
132	135	137	140	142	147	150	158	

Material and Properties:

Material: PP

grammage: 50gsm

Material: Polyester / Nylon / Cotton / Spandex / Velcro

Production: 100% Polyester, 100% Cotton, laboratorial workshop

MOQ: 5000pcs

Order: 5000pcs / 10000pcs / 50000pcs / 100000pcs

Lead time: 80% of the order within 15 days

Payment: T / T, PayPal, L / C, D / A, D / P, Western Union

Service: OEM / ODM

Certification: CE, ISO13485

Capacity: 10,000 / 1 * 40HQ container 10-15 days

Additional information:



orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INTERNATIONAL

