

Product details

Product name: [Placeholder]

Product description: [Placeholder]
Material: [Placeholder]; abrasion-
resistant; [Placeholder]; repellency [Placeholder]

Weight: 25-65gsm

Color: [Placeholder]

1. [Placeholder]
2. [Placeholder]

Lead time: [Placeholder]

1. Sample [Placeholder]: 3-5days
2. [Placeholder]
3. [Placeholder]: 1 [Placeholder]
4. [Placeholder]

Material	PP + PE
Weight	35gsm [Placeholder]
Color	[Placeholder]
Size	120x130cm, 130x150cm, 150x160cm, 150x190cm, 170x170cm
Usage	[Placeholder] / [Placeholder] / [Placeholder]
Quantity	[Placeholder], 10 [Placeholder] / [Placeholder] 10 [Placeholder] / [Placeholder] [Placeholder] 1 [Placeholder] / [Placeholder] 100 [Placeholder] / [Placeholder]
Lead time	[Placeholder]
Dimensions	58X28X48 [Placeholder].

Product display



orient Honest
MEDICAL INTERNATIONAL






orient Honest
MEDICAL INSTRUMENTS



Our business process



The regular shipment way



FAQ

Q: 1. Are your products available for export to the United States?

A: Yes, we have been exporting since 2003. We have a long history of serving customers in the United States since 2013.

Q: 2. What is your minimum order quantity (MOQ)?

A: Our MOQ is USD18 per unit, with a 20% discount for orders over 100 units.

Q: 3. How long does it take to receive the goods?

A: For orders of 380-400 units, it takes 10 QC (Quality Control) days to receive the goods. For orders of 10 Int'l (International) units, it takes 3-8 days to receive the goods. For Int'l orders, it takes 10 days to receive the goods.

Q: 4. What is your payment method?

A: We accept 90% payment in advance and 10% payment upon delivery.