

[KF94](#) [KF94](#) [4](#)

>

KF94 4-layer non-woven face mask with earloop, CE, ISO13485, FDA approved. Features include Super Soft Earloop and 4-layer construction.

- Super Soft Earloop

- 4-layer construction

- Super Soft Earloop

>>

[KF94](#) [KF94](#) [4](#) >

Material	Layer 1: PP Non Woven (Polypropylene) Layer 2: Non-woven fabric Layer 3: Non-woven fabric Layer 4: PP Non-woven fabric
Structure	4ply
Size (IZE)	20*8 CM
Weight (W)	30GSM + 25GSM + 25GSM+25GSM, C OULD Be
Standard	CE, ISO13485, FDA
Moq	20000pcs
Lead time	5-7 days
Quantity	50 PCS/B OX, 2000 S/CTN
Color	White
Price	7-30 USD
Notes	Non-woven fabric, Polypropylene, Super Soft Earloop, 4-layer construction

[KF94](#)



orient Honest
MEDICAL INTERNATIONAL



orient Honest
MEDICAL INTERNATIONAL









